

Noninvasive Cardiovascular Monitoring

CPT Fletcher M Boulware Jr.
MD

Grand Rounds Presentation

Case Presentation

- 80 YOF with CAD and poorly controlled DM admitted to the ICU for sepsis from a left foot cellulitis
- She developed hypoxemia and worsening pulmonary infiltrates and was subsequently intubated
- Aggressive fluid resuscitation along with vasopressor therapy failed to produce adequate perfusion pressures

Case continued

- Decision made to place PAC to guide therapy
- PAC placed without complication, yet 15 minutes after balloon inflation massive hemoptysis noted from ET tube
- Portable CXR shows PAC tip into distal pulmonary artery

07-08-21

E 5-1

F: 508

A

E

V A H C MENT PL 33125

PULH 3 PER 853
24-APR-02
1-13-27
Page 14
PUL- 450
ALL- 157



07-68-11

R 0-1

L 0-1

PULH 3 PER 600

24-APR-68

1:32:04

Page 27


R-0: 175

L-0: 55



E

U.S. AIR FORCE HOSPITAL FL 33185



Noninvasive technology

- Procedures
- Anesthetic techniques / Drugs
- Imaging
- Monitoring

Noninvasive cardiovascular monitoring

- Continuous arterial blood pressure
- Central venous pressure / PAOP
- Cardiac output
- Novel measurements

ABP

- Techniques
 - Auscultation
 - Oscillometry
 - Tonometry

Indications for continuous ABP

- Large fluctuations in BP
- Non-pulsatile flow – CPB
- Frequent lab draws

Contraindications

- Need adequate perfusion
- Infection
- Hypocoagulable state
- Dialysis grafts

Noninvasive ABP



Medwave Vasotrac



Medwave Primo



Physics

$$\text{pressure} = \frac{\text{force}}{\text{area}}$$

Noninvasive ABP

- Contraindications
- Advantages
- Disadvantages

Accuracy - T-Line

- 25 patients undergoing GA
- Compared to invasive A-line in contralateral arm
- The mean \pm sd bias (mm Hg)
 - Sys - 1.7 ± 7.0
 - Diastolic - 2.3 ± 6.9
 - Mean - 1.7 ± 5.3
- Range systolic 41 to 189

Gregory M. Janelle and Nikolaus Gravenstein **An Accuracy Evaluation of the T-Line® Tensymeter (Continuous Noninvasive Blood Pressure Management Device) versus Conventional Invasive Radial Artery Monitoring in Surgical Patients** Anesth. Analg. 2006 102: 484-490

More accuracy data

- Nearly 150,000 paired data points over more than 57 hours of monitoring were collected and analyzed
- Systolic Pressure - Mean Error - 1.0 - STD 7.2
- Diastolic Pressure - Mean Error - 0.8 STD 7.5
- Mean Pressure - Mean Error - 1.3 STD 6.2
- AAMI SP-10:2002 - Accuracy Specification of a mean error not exceeding ± 5 mm Hg and a standard deviation not exceeding 8 mm Hg


Vasotrac Accuracy

- 80 Patients for surgery or in ICU
- Vasotrac vs contralateral radial artery tonometer
- R values
 - Systolic=0.93
 - Diastolic=0.89
 - Mean=0.95
- Pulse Rate 0.95

A New Noninvasive Method to Measure Blood Pressure: Results of a Multicenter Trial.

Anesthesiology. 91(3):686, September 1999. Belani, Kumar M.D. *; Ozaki, Makoto M.D. +; Hynson, James M.D. ++; Hartmann, Thomas M.D. [S]; Reyford, Hugo M.D. [//]; Martino, Jean-Marc M.D. #; Poliac, Marius Ph.D. **; Miller, Ronald M.D. ++

61 YOF for Distal Bypass



Obesity T-line

- 9 Female patients
- BMI > 30
- Systolic 10-beat Average Error – Mean 1.6
STD 7.6
- Diastolic 10-beat Average Error – Mean
4.2 STD 8.1
- Mean 10-beat Average Error - Mean 3.3
STD 7.5

Evaluation of Safety and Accuracy of the T-Line® Tensymeter (Continuous Non-Invasive Blood Pressure Management Device) versus Conventional Invasive Radial Artery Tonometry in Clinically Obese Surgical Patients R. L. Marcus, S. Ahmad, R. Glassenberg, P. Fitzgerald; (Northwestern University Feinberg School of Medicine, Chicago, IL)

Vasotrac - Obesity

- 33 Patients for gastric banding or gastric bypass
- Vasotrac vs A-line vs NIBP Cuff
- VS A-line
 - Sys 0.86
 - Diastolic 0.85
 - Mean 0.89
- VS NIBP Cuff
 - Sys - 0.60
 - Dia - 0.59
 - Mean 0.34

Continual Non-Invasive Blood Pressure Monitoring with the Vasotrac™
-Experience in the Morbidly Obese Authors: DS Beebe MD, I Ostanniy MD, V Komanduri MS, M Poliac PhD, KG Belani MBBS MSPublication/Journal: ANESTHESIOLOGY, V89, No. 3A, September 1998

Summary noninvasive ABP

- Accurate
- Useful?
- Expensive
- When will we see



Cardiac Output

- $CO = HR \times SV$
- Normal in 70kg is about 5L/min
- $Flow = Pressure / Resistance$
- $CO = BP / SVR$

Cardiac Output

- Stewart-Hamilton

$$\dot{Q} = \frac{I}{\int_0^{\infty} C_1 dt} \quad (1)$$

Indicators

- Dyes
- Radioactive tracers
- Thermodilution
- Lithium dilution

Thermodilution

$$\dot{Q} = \frac{(T_B - T_I) \cdot K}{\int_0^{\infty} \Delta T_B(t) dt} \quad (2)$$

Errors

- Shunts
- Right sided valvular lesions
- Inadequate delivery of thermal indicator
- Central venous injection site within the catheter introducer sheath – malposition
- Warming of iced injectate
- Thermistor malfunction from fibrin or clot
- After cardiopulmonary bypass – pulmonary blood temperature change
- Rapid intravenous fluid administration
- Respiratory cycle influences

Fick

$$\dot{Q} = \frac{\dot{V}O_2}{(CaO_2 - C\bar{v}O_2) \cdot 10} \quad (3)$$

$$S\bar{v}O_2 = SaO_2 - \frac{\dot{V}O_2}{\dot{Q} \cdot 1.36 \cdot Hgb} \quad (4)$$

Indications

- To optimize fluid/pharmacologic management in patients with abnormal cardiac or vascular function to achieve normal or supernormal oxygen delivery to organs

Complications

- Pulmonary artery rupture
- Right heart block
- Cardiac puncture
- Dysrhythmias
- Knotting
- Air embolus
- Infection / Sepsis

Less invasive cardiac output

- Doppler
- Partial inert gas rebreathing
- Intravascular pulse contour analysis
 - PiCO
 - PulseCO/LiDCO

Doppler

- Transtracheal
- Esophageal
- TEE

Partial inert gas rebreathing


$$\dot{Q} = \frac{\dot{V}\text{CO}_2}{(\bar{C}\text{VCO}_2 - \text{CaCO}_2)} \quad (14)$$

PiCO

- Uses arterial waveform analysis
- Requires transpulmonary dilution (central line to fem /axillary catheter) for calibration

LiDCO / Pulse CO

- Uses existing arterial line
- Analysis pulse waveform to give a raw measurement of volume changes in the arterial system
- Lithium chloride calibration



Validation

- 19 cardiac surgery patients in ICU
- Pulse contour analysis vs pulmonary thermodilution
- Mean difference 0.31 +/- SD 1.25 L/min
- $R = 0.88$

J Cardiothorac Vasc Anesth. 2000 Apr;14(2):125-9 Beat-to-beat measurement of cardiac output by intravascular pulse contour analysis: a prospective criterion standard study in patients after cardiac surgery. Zollner C, Haller M, Weis M, Morstedt K, Lamm P, Kilger E, Goetz AE.

No good in Off-pump CABG

- 23 Patients for OPCABG
- PulseCO vs thermodilution during
 - Sternotomy – $r^2 = 0.49$
 - Mediastinal opening – $r^2 = 0.52$
 - End of surgery – $r^2 = 0.55$

Can J Anaesth. 2005 May;52(5):530-4 Cardiac output by PulseCO is not interchangeable with thermodilution in patients undergoing OPCAB. Yamashita K, Nishiyama T, Yokoyama T, Abe H, Manabe M


Disadvantages

- Must have functioning A-line
- Calibration required every eight hours
- Not accurate in all clinical scenarios

Noninvasive cardiac output

- Bioimpedance
- Flotrac


Bioimpedance (BioZ)



Bioimpedance

$$SV = \frac{\rho L^2}{Z_0^2} \cdot VET \cdot \max \frac{dZ}{dt} \quad (13)$$

ECG and BioZ waveforms



ICG NORMAL RANGES

05 JUN 18:14

		LOW	HIGH
ICG HR	60	60	100
MAP		70	105
CO	6.6	4.0	8.0
CI	3.3	2.5	4.5
SV	110	60	130
SI	55	35	65
SVR		900	1400
SVRI		1900	2400
TFC	32	30	50
ACI	135	70	150
VI	59	35	65
LVSWI		40	60
LCWI		3.0	5.5
STR	0.36	0.30	0.50
PEP	121		
LVET	334		
eDO2I			

III



V5




160



0



30



0



ICG-Delayed

P

PVC 0

69

ECG
150
50

135/

73

PA2
350
-99

94

22

CVP 3
350
-99

CO

6.6

CI

3.3

ICG

SVR

X

TFC

32

MAIN
MENU

PATIENT
INFO

FAST
LOOK

> NORMAL RANGES
UNITS OF MEASURE
RETURN

TRENDS

HELP

WAVEFORM:
ICG

SPEED:
25

↑

↓

CHECK
LEADS

BioZ Accuracy

- BioZ vs PA Thermodilution after cardiac surgery
- 20 Post CPB patients
- Simultaneous measurements
- Correlation Coefficient of 0.99

J Cardiothorac Vasc Anesth. 2002 Feb;16(1):8-14. Equivalence of bioimpedance and thermodilution in measuring cardiac index after cardiac surgery. **Sageman** WS, **Riffenburgh** RH, **Spiess** BD Pulmonary Medicine, Monterey, CA , USA

BioZ Accuracy

- 53 Post CPB ICU patients
- 210 pairs of measurements
- BioZ vs TD - $R = 0.811$
- Variability in each technique
6.3% vs 24%

BioZ in CHF

- 33 Stable ICU patients with Dx of CHF
- 4 dropped due to inability to get reliable measurements




Figure 1 Concordance correlation plot for cardiac output. Diagonal line shows where data points would be if both techniques provided identical values.

Equivalence of Bioimpedance
and Thermodilution in Measuring Cardiac Output and Index in Patients with Advanced,

BioZ in pulmonary HTN

- 39 Patients referred for eval of PHTN
- Right heart cath
- BioZ vs TD vs Fick
- Correlation
 - BioZ vs TD – 0.80
 - BioZ vs Fick – 0.84

Comparison of Impedance Cardiography to Direct Fick and Thermodilution Cardiac Output Determination in Pulmonary Arterial Hypertension. *Congest Heart Fail.* 2004;10(2 suppl 2):7-10. Authors: Yung GL, Fedullo PF, Kinninger K, Johnson W, Channick RN.

Eliminate need for PAC

- 107 patients in coronary care unit
- 14 patients determined to need hemodynamic monitoring by treating physician
- Physicians given data from BioZ
- 10/14 patients determined not to need a PAC
- Physicians said the BioZ was helpful in all 10
- Physicians said that it improved outcome in 6/10 patients

Advantages

- Easy to initiate
- Limited clinical expertise needed
- Noninvasive
- Accurate
- Outpatient / Emergency
- Cost effective

Disadvantages

- Not accurate in all clinical scenarios
- New / Unfamiliar

Flotrac

- Connects to existing A-line
- No calibration method required



Validation

- 84 Surgical / ICU patients
- Already had PAC in place
- Flotrac vs thermodilution
- Grouped measurements - 562 data points
- Mean difference 0.19 SD 1.28 L/min had

Validation of Continuous Cardiac Output Measurement Using Arterial Pressure Waveforms William T. McGee et al Critical Care Supplement Mar 05

Disadvantages to less/noninvasive cardiac output


- PA pressures
- No PA occlusion pressure
- Mixed venous saturation

Summary noninvasive CO

- Accuracy
- Utility
- Cost


Novel cardiovascular parameters

Hemonix HemoView PeriVasc




Non-invasive intermittent CVP

V A B C T 0 6 - 0 1 0 2 0 4 2 5 . h e m u p



Non-invasive CVP validation

Blinded Hemoview Results (80 Readings on 8 Subjects)
Correlation Coefficient = .82



Intravascular volume

- Ratio of arterial blood to venous blood volume to total blood volume
- Ratio of intravascular blood volume to extravascular blood volume
- Ratio of intravascular blood volume to total body volume?!

Other measured parameters

- CO
- Continuous ABP – without any pressure applied
- Hematocrit
- \$\$\$

And beyond

- Existing technology implemented
- Pharmacogenomics
- TCI
- Nanomedicine – Respirocytes

References

Fegler G: Measurement of cardiac output in anesthetized animals by a thermodilution method. Q J Exp Physiol 39:153-164, 1954.

Berger RL, Weisel RD, Vito L, et al: Cardiac output measurement by thermodilution during cardiac operations. Ann Thorac Surg 21:43-47, 1976.

Jansen JRC: The thermodilution method for the clinical assessment of cardiac output. Intensive Care Med 21:691-697, 1995.

Levett JM, Replogle RL: Thermodilution cardiac output: A critical analysis and review of the literature. J Surg Res 27:392-404, 1979.

Journal of Intensive Care Medicine, Vol. 19, No. 5, 291-296 (2004)
DOI: 10.1177/0885066604265255 Pulmonary Artery Rupture Induced by a Pulmonary Artery Catheter: A Case Report and Review of the Literature Alexandre R. Abreu, MD Michael A. Campos, MD Division of Pulmonary and Critical Care Medicine, Department of Medicine, University of Miami School of Medicine, Miami, Florida Bruce P. Krieger, MD

Nishikawa T, Dohi S: Errors in the measurement of cardiac output by thermodilution. Can J Anaesth 40:142-153, 1993.

Evaluation of Safety and Accuracy of the T-Line® Tensymeter
(Continuous Non-Invasive Blood Pressure Management Device) versus Conventional Invasive Radial Artery
Tonometry in Clinically Obese Surgical Patients R. L. Marcus, S. Ahmad, R. Glassenberg, P. Fitzgerald; (Northwestern University Feinberg School of Medicine, Chicago, IL) Presented at the International Anesthesia Research Society 80th Clinical & Scientific Congress. San Francisco, CA March 24-28, 2006. Presentation Number: S-154

References continued

An Accuracy Evaluation of the T-Line® Tensymeter
(Continuous Noninvasive Blood Pressure Management Device) versus Conventional Invasive Radial Artery Monitoring in Surgical Patients

Gregory M. Janelle, M.D., and Nikolaus Gravenstein, M.D. (Department of Anesthesiology, University of Florida College of Medicine, Gainesville, FL)

Reliability of Hypotension Detection with Noninvasive Radial Artery Beat-to-Beat versus Upper Arm Cuff BP Monitoring

Ron Dueck, M.D. (Anesthesiology, University of California, San Diego and VA San Diego Healthcare System, San Diego, CA), Leslie C. Jameson, M.D. (Anesthesiology, University of Colorado Health Science Center, School of Medicine, Denver, CO) *Presented at the Society for Technology in Anesthesia 2006 Annual Meeting. San Diego, CA January 17-21.*

The New TL-150 Tensymeter

® Continuous Noninvasive BP (CNBP) Versus Direct Radial Artery Measurement of Labile Blood Pressure

Ron Dueck, M.D. (Anesthesiology, University of California, San Diego and VA San Diego Healthcare System, San Diego, CA)

Comparison of NIBP with T-Line® Tensymeter

® (Continuous NIBP Device) for BP Concordance in ASA II-IV Patients Leslie C. Jameson, M.D., Colleen K. Dingmann, R.N., Ph.D. (Anesthesiology, University of Colorado Health Science Center, School of Medicine, Denver, CO)

Temporal Performance of the T-Line® Tensymeter

® (Continuous Non-Invasive Blood Pressure Management Device) versus Conventional Invasive Radial Artery

Tonometry in Surgical Patients. G. M. Janelle, S. Gallant, A. Butler, M. Wilson, N. Gravenstein (University of Florida, Gainesville, FL; Tensys® Medical Inc., San Diego, CA)

References continued

Non-Invasive, Beat-to-Beat Radial Arterial Pressure Monitor (T-Line® TL-100) Provides Blood Pressure Equivalent to Standard Intra-Arterial Catheter

Leslie C. Jameson, M.D. (Anesthesiology, University of Colorado, School of Medicine, Denver, CO)

Comparison of the T-Line® Tensymeter

® (Continuous Non-Invasive Blood Pressure Management Device) with Conventional Invasive Radial Artery

Tonometry in Surgical Patients Gregory M. Janelle, M.D., Anthony Butler, Mark Wilson, Nikolaus Gravenstein, M.D. (Anesthesiology, University of Florida, Gainesville, FL)

Continuous Blood Pressure Monitoring and Patient Safety David B. Swedlow, M.D. (formerly Senior Vice President of Medical Affairs and Technology, Nellcor Puritan Bennett, Inc.) *Data on file, Tensys Medical, Inc*

A Comparison of the T-Line® Tensymeter® Device vs

Direct Radial Arterial Pressure during Major Surgery Thomas Bellehumeur, M.D., James Ramsay, M.D. (Department of Anesthesiology, Emory University School of Medicine, Atlanta, GA)

Accuracy of the T-Line® Tensymeter

® (Continuous Non-Invasive Blood Pressure Management Device) Versus Conventional Invasive Radial Artery

Tonometry in Surgical Cases with Induced Hypotension Peter Szmuk, M.D., Evan Pivalizza, M.D., Ralf Gebhard, M.D., Didier Sciard, M.D., Robert D. Warters, M.D. (Anesthesiology, University of Texas Medical School at Houston and Outcome Research Institute, Houston, TX)

PeriVasc: A Noninvasive Vascular Microscope Patent # 6,749,567 other patents pending Charles L. Davis Hemonix, Inc.

J Cardiothorac Vasc Anesth. 2002 Feb;16(1):8-14. Equivalence of bioimpedance and thermodilution in measuring cardiac index after cardiac surgery. **Sageman WS, Riffenburgh RH, Spiess BD** Pulmonary Medicine, Monterey, CA , USA

References continued

Equivalence of Bioimpedance and Thermodilution in Measuring Cardiac Output and Index in Patients with Advanced, Decompensated Chronic Heart Failure Hospitalized in Critical Care *Am J Crit Care*. 2004;13(6):469-479. Authors: Albert N, Hail M, Li J, Young JB.

Impedance Cardiography: The Next Vital Sign Technology? *Chest*. 2003;123(6):2028-2033. Authors: Van De Water JM, Miller TW.

Comparison of Impedance Cardiography to Direct Fick and Thermodilution Cardiac Output Determination in Pulmonary Arterial Hypertension. *Congest Heart Fail*. 2004;10(2 suppl 2):7-10. Authors: Yung GL, Fedullo PF, Kinninger K, Johnson W, Channick RN.

Congest Heart Fail. 2004 Mar-Apr;10(2 Suppl 2):17-21 Evaluation of impedance cardiography as an alternative to pulmonary artery catheterization in critically ill patients. Silver MA, Cianci P, Brennan S, Longeran-Thomas H, Ahmad F.

Can J Anaesth. 2005 May;52(5):530-4 Cardiac output by PulseCO is not interchangeable with thermodilution in patients undergoing OPCAB. Yamashita K, Nishiyama T, Yokoyama T, Abe H, Manabe M.


J Cardiothorac Vasc Anesth. 2000 Apr;14(2):125-9 Beat-to-beat measurement of cardiac output by intravascular pulse contour analysis: a prospective criterion standard study in patients after cardiac surgery. Zollner C, Haller M, Weis M, Morstedt K, Lamm P, Kilger E, Goetz AE.

Validation of Continuous Cardiac Output Measurement Using Arterial Pressure Waveforms
William T. McGee et al *Critical Care Supplement* Mar 05


Two headed catheter

Medscape®

www.medscape.com



Nice PAC“KNOT”



Balloon leak

